

SKIN REJUVENATION

Vital Injector2, the new level of drug-delivery injection system which will be your cornerstone for the future

Vital Injector2 represents a smart and strategic investment for the clinic so that you can gain a strong competitive differentiator

- Engineered with unparalleled hygiene-centric system - detachable single-use filter
- Praised for accurately adjusted depth control, powered by vacuum assisted technology
- Reduced drug-loss during injection procedure
- Improved operational efficiency

SPECIFICATION

Needle Information	31 G, 5 pin Multi Needle Type
Needle Penetration Control	0~5 mm (Interval 0.2mm)
Syringe Type	1cc / 3cc /5cc
Operation Mode	Intermittent, Continuous
Input Power	AC 100 - 240V, 50/60Hz
Dimension (WxDxH mm)	332 X 315 X 210
Weight	9kg



EunSung Global Corp.

#1102~1105, 53, Gasan digital 2-ro,
Geumcheon-gu, Seoul, 153-706, Korea
Tel: +82-2-514-9713(Rep) Fax: +82-2-514-9716
e-mail: es@esglobal.co.kr / www.esglobal.co.kr



Medical Professionals
Recommended

VITAL INJECTOR2™

Surprisingly Vital

Introducing the All-Renewed
Rejuvenating Injector

A grandeur initiative we all share in search of clean and accurate beauty care for the people



Strong performance with innovative accuracy make the Vital Injector2 a highly desirable equipment for a variety of clinical settings

Advanced capabilities

Now, Featuring with Detachable Hygienic - Filter Technology


Built with advanced hygiene-centric technologies inside, this small but highly capable item effectively protect the system to use it safer and cleaner help you make reliable, confident treatment



Incorporating the latest technology vacuum assisted needle

Newly engineered needle delivers all the essential functionality to meet the clinical needs

- **Pressure release** needle enables to inject all needles beneath skin and subcutaneous tissue - minimize drug loss
- **Controlled depth adjustment** from 0mm to 5mm at 0,2mm interval, delivers customized treatment to accurately target sub dermis where most of collagen and elastin exists
- **31G, 5 Pin multi needle** provides dynamically faster treatment



VITAL INJECTOR2™ is your guide to more professional treatment.


Engineered for operational efficiency



- User-friendly interface simplifies and speeds operation
- 5 inch LCD monitor
- Various sizes of syringes applicable (1 to 5cc) to select from
- Variety of speed formulations allowed
- Convenient audible alarm signal
- activation initiated when drug within the syringe is almost used up

Dermal injection with Vital Injector2 means real dermal hydration

Healthcare professional recommended skin care science that leverages proven formulation at full capacity



- Hydration
- Elasticity
- Roughness
- Smoothness
- Complexion
- Texture

